


Software: Medical Director Version: 3.16b (3.16.16.10323)

Pathology Lab Abnormal Flag is retrieved from the [ABNORMAL_FLAGS] column in the [MD_PATHOLOGY_ATOM] table. The flag displayed is direct from the pathology results, determined by the Pathology Lab. POLAR does not choose the classification. This is populated from the Results tab on the Patient form:



Software: Best Practice

Version: Summit 2014 Build 1.8.4.642

Pathology Lab Abnormal Flag is retrieved from the [ABNORMALFLAG] column in the [REPORTVALUES] table. The flag displayed is direct from the pathology results, determined by the Pathology Lab. POLAR does not choose the classification.

This is populated from Investigation Reports under the Patient form:

